STANDARD[™] LipidoCare - Lipid Test Strip -

Caution

Please follow the instructions in this package insert and the User Manual of the STANDARD[™] LipidoCare Analyzer. Please note that these instructions cover the use of the STANDARD[™] Lipid Profile strips, the STANDARD[™] Total Cholesterol Strips and the STANDARD[™] Triglycerides strips. Please refer to the package to check which specific testing strips you have.

Lipid measurements are used in the diagnosis and treatment of lipoprotein metabolism, lipid disorders, atherosclerosis, and various renal and liver diseases. The determination of cholesterol provides the physician with initial information regarding the state of the patients' lipid metabolism. This information serves as a basis for further diagnostic measures, therapeutic decisions or monitoring the progress of the patient. However, self-monitoring for this system should not be used for diagnosis. The results of self-monitoring should be recorded and discussed with the physician.

INTENDED USE

For the quantitative determination of total cholesterol, High-Density Lipoprotein (HDL) cholesterol and Triglycerides in the blood. Estimated values for Low-Density Lipoprotein (LDL) cholesterol, LDL/HDL ratio and non-HDL cholesterol are calculated by the STANDARD[™] LipidoCare Analyzer and Lipid Profile Test Strip. This system is intended for self-testing as well as professional use.

PRINCIPLES OF THE TEST

STANDARD[™] LipidoCare system combines enzymatic methodology[™] and solid-phase technology to measure Total cholesterol, HDL cholesterol and Triglycerides. Samples used for testing can be whole blood from both capillary (for self-testing and professional use) and venous (for professional use) and serum or plasma. When a blood sample is applied to the STANDARD LipidoCare Lipid Test Strip, the blood reacts to produce a colour that is read by using the reflectance photometry of the analyzer. The level of colour produced is proportional to the concentration. The enzymatic reactions that occur are listed below.

Total cholesterol (TC)

Cholesterol ester $\xrightarrow{\text{cholesterol esterase}}$ Cholesterol + Fatty acid

Cholesterol + O₂ $\xrightarrow{\text{cholesterol oxidase}}$ Cholesterol-3-one + H₂O₂

2H₂O₂ + 4-AAP + Coupler $\xrightarrow{\text{peroxidase}}$ Quinoneimine + 4H₂O
(Toluidine) (Colored-dye)

HDL cholesterol

Total cholesterol specimen $\xrightarrow{\text{VLDL, LDL, HDL depleting reagent}}$ VLDL, LDL depleted specimen

Cholesterol Esters $\xrightarrow{\text{cholesterol esterase}}$ Cholesterol + Fatty acid

Cholesterol + O₂ $\xrightarrow{\text{cholesterol oxidase}}$ Cholesterol-3-one + H₂O₂

2H₂O₂ + 4-AAP + Coupler $\xrightarrow{\text{peroxidase}}$ Quinoneimine + 4H₂O
(Toluidine) (Colored-dye)

Triglycerides (TG)

Triglycerides $\xrightarrow{\text{lipoprotein lipase (LPL)}}$ Glycerol + Fatty Acids

Glycerol + ATP $\xrightarrow{\text{glycerol kinase (GK)}}$ Glycerol-3-phosphate + ADP

Glycerol-3-Phosphate + O₂ $\xrightarrow{\text{glycerol phosphate oxidase (GPO)}}$ 3-dihydroxyacetone-P₂ + H₂O

H₂O₂ + 4-AAP + Coupler $\xrightarrow{\text{peroxidase}}$ Quinoneimine + HCl + 2H₂O
(Toluidine) (Colored-dye)

REAGENT COMPOSITION

Each STANDARD LipidoCare Lipid Test Strip contains the following active ingredients.

Cholesterol Esterase (Microorganism)	≥ 3 Unit
Cholesterol Oxidase (Microorganism)	≥ 1.8 Unit
Peroxidase (Biosensitizer)	≥ 0.75 Unit
4-aminoantipyrine	≥ 37 µg
Toluidine	≥ 50 µg
Dextran Sulfate	≥ 0.16 mg
Glycerol-3-Phosphate Oxidase (Microorganism)	≥ 0.24 Unit
Glycerol Kinase (Microorganism)	≥ 3.44 Unit
ATP (Microorganism)	≥ 153 mg
Lipoprotein Lipase (Microorganism)	≥ 64.8 Unit

Each pouch contains a maximum of 1 g desiccant.

PRECAUTIONS AND WARNINGS

- The STANDARD[™] LipidoCare Lipid Test strip should only be used with the STANDARD[™] LipidoCare Analyzer and the STANDARD[™] LipidoCare Plus Analyzer.
- The STANDARD[™] LipidoCare Lipid Test Strip should not be used beyond the printed expiration date.
- Make sure that the number on the codechip matches the code number printed on the strip pouch.
- The STANDARD[™] LipidoCare Lipid Test System should be tested between 18-32°C (64-90°F).
- A test strip is for single use only. Do not reuse.
- Insert a test strip into the test strip slot with the blood application chamber facing up and toward the analyzer.
- Insert a codechip into the codechip slot with the code number facing up and toward the analyzer.
- Ensure the proper volume for the test strip is used. The volume should be 10µl for single test strips and 35µl for multi-test strips (Lipid Profile).
- Insert a test strip into the analyzer gently until it will go no further.
- Only apply the sample to the blood application area of a test strip.
- Do not ingest.
- Discard the used test strip according to local guidelines.

STORAGE AND HANDLING

- The STANDARD[™] LipidoCare Lipid Test strip should be stored in the sealed pouch.
- Strips should ideally be stored between 2-32°C (36-90°F). Test strips may be stored in a refrigerator between 2-8°C (36-46°F), but must be brought to room temperature before use.
- Use a test strip immediately after removing it from the pouch.
- Keep the codechip either in the analyzer or the test strip package.
- Keep away from heat and direct sunlight.

SAMPLE

Sample Volume

- 10µl for single test strip (TC, TG test strip)
- 35µl for Lipid Profile test strip

Sample Type and Handling

Capillary whole blood (For self-testing and professional use)

- Use fresh capillary whole blood only.
- Avoid excess squeezing of the finger, this may lead to inaccurate results.
- Transfer the blood into the test strip immediately after collecting.

Venous blood (For professional use ONLY)

- Use fresh venous whole blood drawn in heparin or EDTA tube, serum or plasma.
- Venous whole blood should be used within 6 hours of collection.
- Samples should be at room temperature for testing.
- Mix all samples by inverting gently 7-8 times before testing.
- Serum and Plasma are also available for professional use.

CAUTION

- Handle and dispose of all materials coming in contact with blood according to universal precautions and guidelines.
- Do not use a tube with any other additives because it may cause inaccurate results.

TEST PROCEDURE

Materials provided

- STANDARD[™] LipidoCare Lipid Test strip
- The codechip related to the lipid test strip
- STANDARD[™] Ezi Tube+ (35 or 10µl) *(cat. no. 02L510B/10C)
- Package Insert

Additional required materials

- STANDARD[™] LipidoCare Analyzer
- STANDARD[™] LipidoCare Plus Analyzer
- Alcohol swabs to clean puncture site
- Lancets for capillary or venous blood collection supplies (vacuum collection tubes, needles and tube holders)
- STANDARD[™] Micropipette (35 or 10µl) and tips
- Gloves
- SDB Lipid Control Solution
- Bio-hazardous waste containers

Coding

Codechip

Whenever a new package of STANDARD[™] LipidoCare Lipid Test strips is opened, the analyzer must be re-coded. If the code number on the package and the code number displayed by the analyzer do not match, the analyzer cannot perform the measurement accurately. To ensure the accuracy of the STANDARD[™] LipidoCare system, you should change the codechip with the new one. This codechip provides the specific information that the analyzer needs to accurately measure your result.

Code Setting

- Make sure the analyzer is turned off. Remove the old codechip if one is installed.
- Insert a new codechip until it snaps into place.
- Turn the analyzer on. A 3-digit code number appears. This number must match the printed code number on your package or on the test strip pouch.

TESTING

CAUTION

For best results, the subject should fast for 6 hours before the sample is collected.


- Insert the codechip that matches the code number on the test strip package or pouch. And then press the ON/OFF button.

- Hold the test strip by the end with the horizontal raised lines. Insert the test strip into the test strip slot until it will go no further. (When the test strip reaches the correct position, the analyzer will beep.)


- Open the measurement chamber flap. When the flashing blood drop symbol appears on the screen, place the blood sample into the blood application hole of the test strip.

NOTE :

- Use STANDARD Ezi Tube + or STANDARD Micropipette for fingerstick samples.
- Use the STANDARD Micropipette for venous blood sample.
- Close the measurement chamber flap.
- In about three minutes, the result will appear on the screen.
- Remove and discard the used strip. Do not add more blood to a test strip that has been used.


CAUTION: The blood drop must fill the entire blood application hole.

MEASURING RANGE

Test strips can measure numerical results in the following ranges :

- Total cholesterol: 100–450 mg/dL (2.59–11.64 mmol/L)
- HDL cholesterol: 25–95 mg/dL (0.65–2.46 mmol/L)
- Triglycerides: 45–650 mg/dL (0.51–7.34 mmol/L)

If the result is outside the measuring range, the STANDARD LipidoCare Analyzer will display 'Lo' or 'Hi' message.

- 'Lo' below the measuring range (LOW)
- 'Hi' above the measuring range (HIGH)

IMPORTANT:
• If you get a result of 'Lo'/'Hi' or an unexpected result for any test, test again with a new test strip.
• Consult your healthcare professional if the result is still 'Hi' or 'Lo'.

Unit

You can change the unit of measurement in the analyzer settings but you can also convert the units (mg/dL to mmol/L or mmol/L to mg/dL), using the appropriate constant.

	mg/dL to mmol/L divide Result mg/dL by	mmol/L to mg/dL multiply Result mmol/L by
Total Cholesterol	38.664	38.664
HDL	38.664	38.664
Triglycerides	88.54	88.54

EXPECTED VALUES

The National Heart, Lung and Blood Institute issued the Third Report of the National Cholesterol Education Program (NCEP) Expert Panel on Detection, Evaluation, and Treatment of High Blood Cholesterol in Adults (Adult Treatment Panel III) in May 2001.¹³ The ATP III report presented the NCEP's updated clinical guidelines for cholesterol testing and management and described the following classifications for cholesterol and triglyceride testing:

	mg/dL	mmol/L	Classification
Total Cholesterol	<200	<5.18	Desirable
	200–239	5.18–6.19	Borderline high
	≥240	≥6.22	High
HDL	<40	<1.03	Low
	≥60	≥1.55	High
	<150	<1.69	Normal
	150–199	1.69–2.25	Borderline high
	200–499	2.26–5.64	High
	≥500	≥5.65	Very high
	<100	<2.59	Optimal
	100–129	2.59–3.34	Near optimal/above optimal
	130–159	3.36–4.11	Borderline high
	160–189	4.14–4.89	High
	≥190	≥4.91	Very high

The ATP III identified HDL levels below 40mg/dL (1.03 mmol/L) as associated with increased risk of coronary heart disease (CHD) in men and women.¹³ A high HDL level above or equal to 60 mg/dL (1.55 mmol/L) is protective and decreases CHD risk.

HDL

HDL is valid only if Triglyceride level is below 650 mg/dL (7.34 mmol/L).

non-HDL

ATP identifies non-HDL (total cholesterol minus HDL) as a secondary target of therapy in persons with high triglycerides (≥200 mg/dL). The goal for non-HDL in persons with high serum triglycerides can be set at 30 mg/dL higher than that for LDL cholesterol on the premise that a VLDL level ≤30 mg/dL is normal.¹³ non-HDL can be calculated using the equation below.

* **non-HDL (calculated) = Total Cholesterol - HDL**

LDL

LDL can be calculated using the equation below. Calculated LDL is an estimation of LDL and valid only if Triglyceride level is 400mg/dL (4.52mmol/L) or below.

Otherwise, the analyzer displays "N/A" for LDL and LDL/HDL.

* **LDL (calculated) = Total Cholesterol - HDL - (Triglycerides/5) (mg/dL)**

* **LDL (calculated) = Total Cholesterol - HDL - (Triglycerides/2.17) (mmol/L)**

QUALITY CONTROL

Quality Control tests are used to ensure that the total system (analyzer, strips, and codechip) is working properly. Users should run a control solution test when results are questionable or to comply with their own facility's quality control requirements. See STANDARD[™] LipidoCare Lipid Test System User Manual for instructions on how to run a control solution test.

LIMITATIONS

- Preservatives: EDTA and heparin in venous blood collection tubes had no effect on the results of the test strip.
- Hematocrits of 30% to 55% for TC and TG test strip and 30 to 52% for Lipid Profile test strip do not affect results.
- Performance of STANDARD LipidoCare Lipid Test System has not been tested on samples from newborns.
- Cosmetics such as handcreams or lotions often contain glycerol. Use of these products may cause inaccurate results.
- Interferences: The following substances may cause false results at levels of the below.

Substances concentration(mg/dL)			
Hemoglobin	> 300	L-Dopa	> 1.6
Ascorbic acid	> 2.5	Urea	> 700
Fructose	> 30	Creatinine	> 30
Gentisic acid	> 1	Glutathione	> 2
Oxytetracycline	> 10	Lactose	TC, TG > 100 / HDL > 50
Cysteine	> 2.5	Bilirubin	TC, HDL > 20 / TG > 10
o-methylodopa	> 2	Nicotinic acid	TC, TG > 30 / HDL > 10
Dipyrrone	> 30	Dopamin	> 2

PERFORMANCE CHARACTERISTICS

1. Accuracy (Method Comparison)

A clinical study was performed by professionals who measured lipid levels on fresh venous whole blood specimens. The results below show that STANDARD[™] LipidoCare Lipid Test strips compared well to a reference method.

Single test strip	No.	Slope	Y-intercept	R
TC	120	0.970	5.099	0.986
TG	120	0.971	4.111	0.995

Lipid Profile test strip	No.	Slope	Y-intercept	R
TC	120	0.994	1.208	0.985
TG	120	0.980	1.305	0.993
HDL	80	0.965	1.662	0.979

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