PRECAUTIONS

- 1) Carefully read the instructions before performing the test. The test is only interpretable if the instructions are carefully respected.
- 2) Keep out of the reach of children.
- 3) Do not use after the expiry date printed on the label and on the protective pouch or if the pouch is damaged.
- 4) Follow strictly the indicated time, blood and diluent quantities.
- 5) Store between +4 °C and +30 °C. Do not freeze.
- 6) Do not re-use the cassette and lancet.
- 7) External use only. DO NOT SWALLOW.
- 8) This test is exclusively intended for use by individuals taking medicines that thin the blood (anticoagulants) or by the individuals with disorders such as hemophilia.
- 10) After use, dispose of all components in accordance with local regulations, seek advice from your pharmacist.

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
CONTENT

THE BOX CONTAINS ALL THE NECESSARY TO PERFORM 1 C-REACTIVE PROTEIN TEST:

- 1 sealed aluminium pouch containing:
 - 1 plastic pipette
 - 1 desiccant bag
- 2 sterile lancets for blood sampling
- 1 disinfectant pad
- 1 instructions leaflet

* Not provided necessary material: absorbent cotton.
 * Only open the protective pouch when you are ready to use the test.
 * The desiccant bag should not be used.

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
The C-Reactive Protein (CRP) is a non specific marker which is mainly produced by the liver and used to diagnose bacterial infectious disease and inflammatory disorders. The CRP is a very sensitive and fast appearing indicator which could therefore be helpful for deciding an antibiotic treatment. Depending on the CRP concentration, different lines will appear allowing the semi-quantitative interpretation of the test result.

In healthy patients, CRP concentration is lower than 8mg/L while the concentration level can go higher than 100mg/L in case of severe infection or during inflammatory process.

Intermediate levels, within 8 and 100mg/L, are concomitant with more or less mildly viral or bacterial infections that can be easily overcome by appropriate treatment ordered by your doctor.

The C-Reactive Protein (CRP) is a non specific marker which is mainly produced by the liver and used to diagnose bacterial infectious disease and inflammatory disorders. The CRP is a very sensitive and fast appearing indicator which could therefore be helpful for deciding an antibiotic treatment. Depending on the CRP concentration, different lines will appear allowing the semi-quantitative interpretation of the test result.

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Cut-off: <8 / 8 - 40 / 40 - 100 / >100 mg/L

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C-REACTIVE PROTEIN TEST
TEST PROTÉINA C-REACTIVA
TEST PROTÉINA C-REATIVA

Self-Test for the detection of C-Reactive Protein levels in blood

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
Legal manufacturer
 Fabricante
 Fabricante

PROCEDURE

- 1) Wash your hands with soap and rinse with clear warm water.
- 2) Tear the protective pouch (from the notch) and only get out the cassette and the pipette. Dispose of the small desiccant bag.
- 3) Unscrew the dropper vial containing the diluent.
- 4) Hold the lancet and carefully rotate 360° the protective cap without pulling it. - **FIG. A**
- 5) Extract and dispose the released cap. - **FIG. B**
- 6) Press the lancet, on the side from where the cap was extracted, against the fingertip (it is advisable the ring finger side). - **FIG. C** The tip retracts automatically and safely after use.
- 7) Keeping the hand down massage the end that was pricked to obtain a blood drop. - **FIG. D**
- 8) **Without pressing the pipette bulb**, put it in contact with the drop. The blood migrates into the pipette through the capillarity to the line indicated on the pipette. - **FIG. E**
- You may massage again your finger to obtain more blood if the line is not reached. As far as possible, avoid of air bubbles.
- 9) Put the blood collected with the pipette into the opened dropper vial, by pressing on the pipette bulb. Press the bulb of the pipette 2 or 3 times to make sure all the blood sample is added into the diluent solution. Then replace the screw cap on the dropper vial and mix well. - **FIG. F**
- 10) **Break the tip of the dropper bottle** containing the blood sample - **FIG. G** and add 4 drops of the diluted sample into the device sample well. - **FIG. H**
- 11) Read the result after 5 minutes.

RESULT INTERPRETATION

**READ THE RESULT AFTER 5 MINUTES.
DO NOT READ THE RESULT AFTER 10 MINUTES.**

The results must be interpreted depending on the intensity colour of the 3 lines that could appear in the device window.

CRP LEVEL LOWER THAN 8 MG/L

Only two lines (no line under the mark T (Test)) appear. This result means that there is neither infection nor inflammation.

CRP LEVEL IN THE RANGE OF 8 TO 40MG/L

Three coloured lines appear in the window but the intensity of the line under T mark is weaker than the intensity of the two other lines. This result means that a viral infection could be on the way. **You should consult your doctor.**


CRP LEVEL IN THE RANGE OF 40 TO 100MG/L

Three coloured lines appear in the window and the intensity of the 3 lines are more or less equal. This result means that a viral or bacterial infection is on the way. **You must consult your doctor.**

CRP LEVEL OVER 100MG/L

Three coloured lines appear in the window but the colour intensity of the lines under T and C marks are stronger than the colour intensity of the middle line. **In case of very high CRP concentration, the middle line may even disappear.** This result means that a severe bacterial infection is on the way. **You must consult your doctor urgently.**

NON VALID RESULT

If no line appears under the mark C (Control), it is not possible to interpret the test, which must be considered as non valid even a coloured line is visible under T (Test) mark and/or the middle line appears. In this case, it is recommended to repeat the test with a new C-REACTIVE PROTEIN Test.

F.A.Q. - QUESTIONS & ANSWERS

HOW DOES C-REACTIVE PROTEIN TEST WORK?

If present in the sample, CRP will interact with the anti CRP antibodies as well as with the CRP antigen placed in the test producing different lines allowing a semi-quantitative measurement of the CRP concentration in the blood sample.

WHEN CAN THIS TEST BE PERFORMED?

C-REACTIVE PROTEIN Test should be performed in the case of acute infection symptoms such as feverishness, fever, headaches or weakness. The test can be performed at any time, but not in subjects younger than 6 years.

CAN THE RESULT BE INCORRECT?

The results are accurate as far as the instructions are carefully respected. Nevertheless, the result can be incorrect if C-REACTIVE PROTEIN Test gets wet before test performing or if the quantity of diluted blood dispensed in the sample well is not sufficient. The plastic pipette provided in the box allows making sure the collected blood volume is correct and is correctly diluted into the diluent dropper vial.

HOW TO INTERPRET THE TEST IF THE COLOUR AND THE INTENSITY OF THE LINES ARE DIFFERENT?

The colour intensity of the lines is important to interpret the test result. So please refer to the pictures and instructions in paragraph "Result interpretation".

WHAT IS THE LINE THAT APPEARS UNDER THE MARK C (CONTROL) FOR?

When this line appears, it only means that the test was per-

forming well.

IF I READ THE RESULT AFTER 10 MINUTES, WILL THE RESULT BE RELIABLE?

No. The result should be read within 10 minutes. The result is reliable up to 10 minutes but should be better read at 5 minutes after adding the sample.

WHAT DO I HAVE TO DO IF THE RESULT IS POSITIVE?


If the result is positive when the CRP is higher than 8mg/L, please refer to the result interpretation and inform your doctor as soon as possible as an antibiotic treatment might be urgently needed.

WHAT DO I HAVE TO DO IF THE RESULT IS NEGATIVE?

If the result is negative when the CRP is lower than 8mg/L, it means that there is neither bacterial nor viral infection on the way. However if the symptoms persist, it is recommended to consult your doctor.

WHAT IS THE ACCURACY OF C-REACTIVE PROTEIN TEST?

The C-REACTIVE PROTEIN Test is accurate and is used for more than 10 years by professionals in the field (hospitals, laboratories...). Evaluation reports show that an overall agreement of at least 98% with reference method. Although this test is reliable, false positive or false negative results can be obtained.

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