

REF OTS-402H	English
--------------	---------

A rapid test for the qualitative detection of human Thyroid Stimulating Hormone (TSH) in whole blood. For self-testing in vitro diagnostic use only.

【SUMMARY】

Thyroid stimulating hormone (also known as thyrotropin, thyrotropic hormone, TSH, or hTSH for human TSH) is a pituitary hormone that stimulates the thyroid gland to produce thyroxine (T₄), and then triiodothyronine (T₃) which stimulates the metabolism of almost every tissue in the body.¹ It is a glycoprotein hormone synthesized and secreted by thyrotrope cells in the anterior pituitary gland, which regulates the endocrine function of the thyroid.^{2,3} TSH (with a half-life of about an hour) stimulates the thyroid gland to secrete the hormone thyroxine (T₄), which has only a slight effect on metabolism. T₄ is converted to triiodothyronine (T₃), which is the active hormone that stimulates metabolism. About 80% of this conversion is in the liver and other organs, and 20% in the thyroid itself.¹

Testing of thyroid stimulating hormone levels in the blood is considered the best initial test for hypothyroidism.⁴ It is important to note the statement from the Subclinical Thyroid Disease Consensus Panel: "There is no single level of TSH at which clinical action is always either indicated or contraindicated. The higher the TSH, the more compelling is the rationale for treatment. It is important to consider the individual clinical context (e.g. pregnancy, lipid profile, ATPO antibodies)."⁵

The TSH Rapid Test Cassette (Whole Blood) is a rapid test that qualitatively detects the presence of TSH in whole Blood specimen at the sensitivity of 5µIU/ml. The TSH Rapid Test Cassette (Whole Blood) is a simple test that utilizes a combination of monoclonal antibodies to selectively detect elevated levels of TSH in whole blood.

【PRECAUTIONS】

Please read all the information in this package insert before performing the test.

- For self-testing in vitro diagnostic use only.
- Do not eat, drink or smoke in the area where the specimens or kits are handled.
- Store in a dry place at 2-30°C (36-86°F), avoiding areas of excess moisture. If the foil packaging is damaged or has been opened, please do not use.
- This test kit is intended to be used as a preliminary test only and repeatedly abnormal results should be discussed with a doctor or medical professional.
- Follow the indicated time strictly.
- Use the test only once. Do not dismantle and touch the test window of the test cassette.
- The kit must not be frozen or used after the expiration date printed on the package.
- Keep out of the reach of children.
- The used test should be discarded according to local regulations.

【STORAGE AND STABILITY】

Store as packaged in the sealed pouch at room temperature or refrigerated (2-30°C). The test is stable through the expiration date printed on the sealed pouch. The test must remain in the sealed pouch until use. **DO NOT FREEZE.** Do not use beyond the expiration date.

【MATERIALS PROVIDED】


- Test cassette
- Capillary Dropper
- Buffer
- Alcohol pad
- Lancets
- Package insert

【MATERIALS REQUIRED BUT NOT PROVIDED】


- Timer

【PROCEDURE】

1. Wash your hands with soap and rinse with clear warm water.
2. Bring the pouch to room temperature before opening it. Open the foil pouch and get out the cassette.
3. Carefully pull off and dispose of the clear plastic cap of the lancet. The lancet is a single use lancet so make sure you do not fire the lancet prematurely.
4. Use the provided alcohol pad to clean the fingertip of the middle or ring finger as the puncture site.
5. Press the lancet (on the side from where the cap was removed) firmly against the fingertip (the side of the ring or middle finger is advised). The tip retracts automatically and safely after use.
6. Keeping the hand down, massage the finger that was pricked to obtain a blood drop.
7. Without squeezing the capillary dropper bulb, gently put it in contact with the edge of the blood drop. The blood should automatically be drawn up the dropper by capillary action. You should fill the dropper up to the black line.
8. You may continue to massage your finger to obtain more blood if the line is not reached. As far as possible, try to avoid air bubbles.
9. Transfer the collected blood into the sample well (S) of the cassette, by squeezing the dropper bulb.
10. Wait for the blood to be totally dispensed in the well. Unscrew the cap of the buffer bottle and add **2 drops of buffer** into the sample well of the cassette.
11. Wait for the coloured line(s) to appear. Read results at **10 minutes**. Results obtained after this time may be inaccurate.




【READING THE RESULTS】



POSITIVE: Two lines appear. Both T (Test) and C (Control) lines appear.
This result means that the TSH level is **higher than the normal (5µIU/mL)** and that you should consult a physician.



NEGATIVE: One line appears. Only control line appears (C).
This result means that the TSH level is not in the range to consider Hypothyroidism.



INVALID: Control line fails to appear. Insufficient specimen volume or incorrect procedural techniques are the most likely reasons for control line failure. Review the procedure and repeat the test with a new test. If the problem persists, discontinue using the test kit immediately and contact your local distributor.

【LIMITATIONS】

- The TSH Rapid Test Cassette (Whole Blood) is for in vitro diagnostic use only. The test should be used for the detection of TSH in whole blood specimens only. Neither the quantitative value nor the rate of increase in TSH concentration can be determined by this qualitative test.
- The TSH Rapid Test Cassette (Whole Blood) is only for screening the primary hypothyroidism of the adult population, not for neonates.
- As with all diagnostic tests, all results must be interpreted together with other clinical information available to the physician.
- A positive test must be confirmed using a quantitative laboratory TSH assay.
- False positive results can occur due to heterophilic (unusual) antibodies. In certain clinical conditions such as central hypothyroidism, TSH levels may be normal/ low, despite hypothyroidism. Medical consultation is recommended to exclude such cases.
- For Central/ Secondary Hypothyroidism, TSH is not a reliable biomarker, which occurs in 1 out of 1,000 Hypothyroidism cases.

【FAQs】

1. How does the TSH test work?

The Thyroid Stimulating Hormone (TSH) activates the thyroid gland. A TSH level over 5µIU/mL (in case of a positive result), indicates an under active thyroid (hypothyroidism), which requires more TSH.

2. When should the test be used?

In cases of hypothyroidism, symptoms such as regularly feeling tired, depressed or cold, weight gain, dry skin, brittle hair, enduring constipation or menstrual cycle irregularities in women, occur and it is recommended to perform a TSH Rapid Test for screening purpose. The TSH Rapid Test can be used any time of the day. However, it cannot and should not be performed in case of hormonal thyroid medical treatment.

3. Can the result be incorrect?

The results are accurate as far as the instructions are carefully respected. Nevertheless, the result can be incorrect if the TSH Rapid Test cassette gets wet before test performing or if the quantity of blood dispensed in the sample well is not sufficient, or if the number of buffer drops are fewer than 2 or more than 3. The capillary dropper provided in the box allows you to make sure the collected blood volume is correct. In addition, due to immunological principles involved, there is the possibility of false results in rare cases. A consultation with a doctor is always recommended for such tests based on immunological principles.

4. How to interpret the test if the colour and the intensity of the lines are different?

The colour and intensity of the lines have no importance for result interpretation. The lines should only be homogeneous and clearly visible. The test should be considered as positive whatever the colour intensity of the test line.

5. If I read the result after 20 minutes, will the result be reliable?

No. The result should be read **10 minutes** after adding the buffer. The result is unreliable after 20 minutes.

6. What do I have to do if the result is positive?

If the result is positive, it means that the TSH level in blood is higher than the normal (5µIU/mL) and that you should consult a physician to show the test result. Then, the physician will decide whether additional analysis should be performed.

7. What do I have to do if the result is negative?

If the result is negative, it means that the TSH level is below 5µIU/mL and is within the normal range. A case of Hyperthyroidism, though rare, cannot be excluded based on such test results. If the symptoms persist, it is recommended to consult a physician.

【BIBLIOGRAPHY】

- Merck Manual of Diagnosis and Therapy, Thyroid gland disorders.
- The American Heritage Dictionary of the English Language, Fourth Edition. Houghton Mifflin Company. 2006. ISBN 0-395-82517-2.
- Sacher R, Richard A. McPherson (2000). Widmann's Clinical Interpretation of Laboratory Tests, 11th ed. F.A. Davis Company. ISBN 0-8036-0270-7.
- So, M; MacIsaac, RJ; Grossmann M (August 2012). "Hypothyroidism". Australian Family Physician 41 (8): 556-62.
- Surks et al. JAMA 291:228, 2004.
- Daniel, GH, Martin, JB, Neuroendocrine Regulation and Diseases of the Anterior Pituitary and Hypothalamus in Wilson, JD, Braunwald, E., Isselbacher, KJ, et. al., Harrison's Principles of Internal Medicine, 12th Edition, McGraw-Hill, Inc., New York, NY, 1991, p. 1666)

Index of Symbols

	Attention, see instructions for use
	For in vitro diagnostic use only
	Store between 2-30°C
	Do not use if package is damaged

	Tests per kit
	Use by
	Lot Number
	Consult Instructions For Use

	Authorized Representative
	Do not reuse
	Catalog #
	Manufacturer

Hangzhou AllTest Biotech Co., Ltd.
#550, Yinhai Street
Hangzhou Economic & Technological Development Area
Hangzhou - 310018, P. R. China
www.alltests.com.cn



EC REP
MedNet GmbH
Borkstrasse 10
48163 Muenster
Germany

Number: 145792600
Effective date: 2017-12-25