Wellion® LEONARDO

Wellion® LEONARDO CHOL Teststreifen

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Wellion® LEONARDO Cholesterol Test trip

IMPORTANT: PLEASE READ THIS INFORMATION AND YOUR WELLION LEONARDO GLU/CHOL USER'S GUIDE BEFORE USING WELLION LEONARDO CHOLESTEROL TEST STRIP.
FOR QUESTIONS AND ASSISTANCE, PLEASE CONTACT THE AUTHORIZED REPRESENTATIVE FOUND AT THE END OF THIS INSERT.

INTENDED USE

The Wellion LEONARDO Cholesterol Test Strip is intended to be used with the Wellion LEONARDO GLU/CHOL Cholesterol Meter to quantitatively measure cholesterol in capillary whole blood drawn from fingertips. The Wellion LEONARDO Cholesterol Test Strip is plasma-calibrated for easy comparison to lab results. It is intended for self-testing by individuals at home or for clinical settings by healthcare professionals.

INTRODUCTION

The Wellion LEONARDO Cholesterol Test Strip measures total cholesterol in blood via enzymatic reaction which produces a current. Blood from a finger tip is drawn into the reaction area which contains the enzyme through passive capillary action. The current produced by the reaction is measured by the meter and is in direct correlation to the concentration of the total cholesterol in the blood sample. The meter shows the test result in 90 seconds.

Reagent Composition

Each cm² of test strip contains:
Cholesterol Esterase 10%
Cholesterol Oxidase 12%
Enzyme Stabilizer 59%
Enzyme Activator 19%

Warnings and Precautions

The Wellion LEONARDO Cholesterol Test Strip is for use outside the body (IN VITRO diagnostic use).

- [1] Do not use test strip after their expiration date. Please check the expiration date on the test strip box. [2] Do not use strips that are wet or damaged. [3] Do not reuse the strips. [4] The test strip should not be used as the basis for diagnosis of disorders involving excess cholesterol in the blood and lipid protein metabolism disorders. [5] Never make any changes in your treatment program or ignore symptoms without consulting your physician. [6] Operation temperature for meter and test strip is between 18°C–38°C (64°F–100°F). [7] You must insert a code chip into the meter for each new box of test strips for proper calibration. The code chip is packaged inside the test strip box and is lot specific. Make sure that both code numbers as shown on the meter after the code chip has been inserted and the code number printed on the test strip box match each other.

Storage and Handling

- [1] Store the test strips between 4°C–30°C (39°F–86°F). Do not freeze. [2] Store out of direct sunlight. [3] Unopened test strips are stable until the expiration date printed on the box when stored properly. [4] The month printed on the label refers to the END of that month. [5] Do not handle the test strip with wet or dirty hands.

Sample Collection and Preparation

Testing must be performed immediately after the sample is obtained.

Test Procedure

See „Testing Your Blood Cholesterol“ in the Wellion LEONARDO GLU/CHOL User's Guide.

Expected Values

Consult with your physician or healthcare professional to determine an appropriate cholesterol target range for you.

Total cholesterol levels:
< 200 mg/dL Desirable
200–239 mg/dL Borderline high
≥ 240 mg/dL High

Reference: 1 "High Blood Cholesterol - What You Need to Know". NHLBI. Publication No. 05-3290

Quality Control

Run Level 1 or Level 2 control solution test. Follow the User's Guide instructions. Do control tests:

- [1] At least once per week to make sure the meter and test strip are working properly. [2] If your test strips were stored at temperature and humidity outside proper storage conditions. [3] When you use your meter for the first time. [4] To practice your testing technique. [5] If you drop your meter.

Use only Wellion LEONARDO Cholesterol Control Solution. Your test results should fall within the control range printed on the test strip bottle. Repeat control solution testing if results fall out of range. Results may fall out of range due to:

- [1] Errors in control solution testing [2] Expired or contaminated control solution [3] Damaged test strip [4] Meter malfunction

IMPORTANT: IF THE RESULT CONTINUES TO FALL OUTSIDE OF THE PRINTED RANGE; CONTACT THE AUTHORIZED REPRESENTATIVE FOUND AT THE END OF THIS INSERT.

Problem Solving

- [1] Confirm the test strip is not expired. [2] Make sure the blood fits the reaction zone. "E 2" will show if there is too little blood. DO NOT add a second drop of blood. Discard the test strip from the meter and retest with a new test strip. [3] Check the system with a control solution test. [4] Refer to "Solving Problems" in the User's Guide for more hints.

Limitations

[1] DO NOT use plasma or serum samples. [2] DO NOT test on neonatal blood samples. [3] DO NOT test on venous or arterial blood samples. [4] Altitudes up to 3150m (10335 feet) will not affect test results. [5] Hematocrit range: 30–55%. [6] High level of ascorbic acid (>1.5 mg/dL) may cause incorrect results. To get correct result, test your cholesterol level in the morning after fasting for at least 8 hours. [7] L-DOPA (>0.5 mg/dL) can give falsely reduced results in therapeutic concentration (1 mg/dL).

The following intererents may affect test results if concentration is greater than therapeutic or physiological concentration:

Interferents	Interference Concentration	Therapeutic Concentration/ Reference Interval
Dopamine	> 0,5 mg/dL	0,03 mg/dL
Methyl-Dopa	> 0,75 mg/dL	0,1–0,75 mg/dL
Triglyceride	> 308 mg/dL	30–300 mg/dL

Performance Evaluation Data

Accuracy

Three lots of Wellion LEONARDO Blood Cholesterol Test Strips were used in the accuracy test of the Wellion LEONARDO GLU/CHOL System.

A capillary blood comparison study between the Wellion LEONARDO CHOL System and Hitachi 717 Chemistry Analyzer yields the following data:

Number of Readings: 327
Sample Range: 101 to 390 mg/dL

Total within +/- 20 %			
314/327 (96,0%)			

Precision

1 lot of Wellion LEONARDO Cholesterol Test StripS was used for within-run repeatability study. Venous blood in heparin-tubes was spiked to 3 kinds of concentration. Blood cholesterol readings were recorded for 1 day resulting in 100 data points for each concentration; as shown in the following tables:

Number of Readings	100	100	100
Mean (mg/dL)	137	193	255
CV%	6,5	5,8	6,3
mCV%	NA	6,2	NA

Wellion® LEONARDO Cholesterol Test trip

IMPORTANT: PLEASE READ THIS INFORMATION AND YOUR WELLION LEONARDO GLU/CHOL USER'S GUIDE BEFORE USING WELLION LEONARDO CHOLESTEROL TEST STRIP.
FOR QUESTIONS AND ASSISTANCE, PLEASE CONTACT THE AUTHORIZED REPRESENTATIVE FOUND AT THE END OF THIS INSERT.

INTENDED USE

The Wellion LEONARDO Cholesterol Test Strip is intended to be used with the Wellion LEONARDO GLU/CHOL Cholesterol Meter to quantitatively measure cholesterol in capillary whole blood drawn from fingertips. The Wellion LEONARDO Cholesterol Test Strip is plasma-calibrated for easy comparison to lab results. It is intended for self-testing by individuals at home or for clinical settings by healthcare professionals.

INTRODUCTION

The Wellion LEONARDO Cholesterol Test Strip measures total cholesterol in blood via enzymatic reaction which produces a current. Blood from a finger tip is drawn into the reaction area which contains the enzyme through passive capillary action. The current produced by the reaction is measured by the meter and is in direct correlation to the concentration of the total cholesterol in the blood sample. The meter shows the test result in 90 seconds.

Reagent Composition

Each cm² of test strip contains:
Cholesterol Esterase 10%
Cholesterol Oxidase 12%
Enzyme Stabilizer 59%
Enzyme Activator 19%

Warnings and Precautions

The Wellion LEONARDO Cholesterol Test Strip is for use outside the body (IN VITRO diagnostic use).

- [1] Do not use test strip after their expiration date. Please check the expiration date on the test strip box. [2] Do not use strips that are wet or damaged. [3] Do not reuse the strips. [4] The test strip should

not be used as the basis for diagnosis of disorders involving excess cholesterol in the blood and lipid protein metabolism disorders. [5] Never make any changes in your treatment program or ignore symptoms without consulting your physician. [6] Operation temperature for meter and test strip is between 18°C–38°C (64°F–100°F). [7] You must insert a code chip into the meter for each new box of test strips for proper calibration. The code chip is packaged inside the test strip box and is lot specific. Make sure that both code numbers as shown on the meter after the code chip has been inserted and the code number printed on the test strip box match each other.

Storage and Handling

- [1] Store the test strips between 4°C–30°C (39°F–86°F). Do not freeze. [2] Store out of direct sunlight. [3] Unopened test strips are stable until the expiration date printed on the box when stored properly. [4] The month printed on the label refers to the END of that month. [5] Do not handle the test strip with wet or dirty hands.
- </

Wellion® LEONARDO Cholesterol Test strip

IMPORTANT: PLEASE READ THIS INFORMATION AND YOUR WELLION LEONARDO GLU/CHOL USER'S GUIDE BEFORE USING WELLION LEONARDO CHOLESTEROL TEST STRIP. FOR QUESTIONS AND ASSISTANCE, PLEASE CONTACT THE AUTHORIZED REPRESENTATIVE FOUND AT THE END OF THIS INSERT.

INTENDED USE

The Wellion LEONARDO Cholesterol Test Strip is intended to be used with the Wellion LEONARDO GLU/CHOL Cholesterol Meter to quantitatively measure cholesterol in capillary whole blood drawn from fingertips. The Wellion LEONARDO Cholesterol Test Strip is plasma-calibrated for easy comparison to lab results. It is intended for self-testing by individuals at home or for clinical settings by healthcare professionals.

INTRODUCTION

The Wellion LEONARDO Cholesterol Test Strip measures total cholesterol in blood via enzymatic reaction which produces a current. Blood from a finger tip is drawn into the reaction area which contains the enzyme through passive capillary action. The current produced by the reaction is measured by the meter and is in direct correlation to the concentration of the total cholesterol in the blood sample. The meter shows the test result in 90 seconds.

Reagent Composition

Each cm² of test strip contains:
 Cholesterol Esterase 10%
 Cholesterol Oxidase 12%
 Enzyme Stabilizer 59%
 Enzyme Activator 19%

Warnings and Precautions

The Wellion LEONARDO Cholesterol Test Strip is for use outside the body (IN VITRO diagnostic use).

[1] Do not use test strip after their expiration date. Please check the expiration date on the test strip box. [2] Do not use strips that are wet or damaged. [3] Do not reuse the strips. [4] The test strip should not be used as the basis for diagnosis of disorders involving excess cholesterol in the blood and lipid protein metabolism disorders. [5] Never make any changes in your treatment program or ignore symptoms without consulting your physician. [6] Operation temperature for meter and test strip is between 18°C-38°C (64°F-100°F). [7] You must insert a code chip into the meter for each new box of test strips for proper calibration. The code chip is packaged inside the test strip box and is lot specific. Make sure that both code numbers as shown on the meter after the code chip has been inserted and the code number printed on the test strip box match each other.

Storage and Handling

[1] Store the test strips between 4°C-30°C (39°F-86°F). Do not freeze. [2] Store out of direct sunlight. [3] Unopened test strips are stable until the expiration date printed on the box when stored properly. [4] The month printed on the label refers to the END of that month. [5] Do not handle the test strip with wet or dirty hands.

Sample Collection and Preparation

Testing must be performed immediately after the sample is obtained.

Test Procedure

See „Testing Your Blood Cholesterol“ in the Wellion LEONARDO GLU/CHOL User's Guide.

Expected Values

Consult with your physician or healthcare professional to determine an appropriate cholesterol target range for you.

Total cholesterol levels:

< 200 mg/dL	Desirable
200-239 mg/dL	Borderline high
≥ 240 mg/dL	High

Reference: 1 "High Blood Cholesterol - What You Need to Know". NHLBI. Publication No. 05-3290

Quality Control

Run Level 1 or Level 2 control solution test. Follow the User's Guide instructions. Do control tests:

[1] At least once per week to make sure the meter and test strip are working properly. [2] If your test strips were stored at temperature and humidity outside proper storage conditions. [3] When you use your meter for the first time. [4] To practice your testing technique. [5] If you drop your meter.

Use only Wellion LEONARDO Cholesterol Control Solution. Your test results should fall within the control range printed on the test strip bottle. Repeat control solution testing if results fall out of range. Results may fall out of range due to:

[1] Errors in control solution testing [2] Expired or contaminated control solution [3] Damaged test strip [4] Meter malfunction

IMPORTANT: IF THE RESULT CONTINUES TO FALL OUTSIDE OF THE PRINTED RANGE; CONTACT THE AUTHORIZED REPRESENTATIVE FOUND AT THE END OF THIS INSERT.

Problem Solving

[1] Confirm the test strip is not expired. [2] Make sure the blood fills the reaction zone. "E 2" will show if there is too little blood. DO NOT add a second drop of blood. Discard the test strip from the meter and retest with a new test strip. [3] Check the system with a control solution test. [4] Refer to "Solving Problems" in the User's Guide for more hints.

Limitations

[1] DO NOT use plasma or serum samples. [2] DO NOT test on neonatal blood samples. [3] DO NOT test on venous or arterial blood samples. [4] Altitudes up to 3150m (10335 feet) will not affect test results. [5] Hematocrit range: 30-55%. [6] High level of ascorbic acid (>1,5 mg/dL) may cause incorrect results. To get correct result, test your cholesterol level in the morning after fasting for at least 8 hours. [7] L-DOPA (>0,5 mg/dL) can give falsely reduced results in therapeutic concentration (1 mg/dL).

The following interferents may affect test results if concentration is greater than therapeutic or physiological concentration:

Interferents	Interference Concentration	Therapeutic Concentration/ Reference Interval
Dopamine	> 0,5 mg/dL	0,03 mg/dL
Methyl-Dopa	> 0,75 mg/dL	0,1-0,75 mg/dL
Triglyceride	> 308 mg/dL	30-300 mg/dL

Performance Evaluation Data

Accuracy

Three lots of Wellion LEONARDO Blood Cholesterol Test Strips were used in the accuracy test of the Wellion LEONARDO GLU/CHOL System. A capillary blood comparison study between the Wellion LEONARDO CHOL System and Hitachi 717 Chemistry Analyzer yields the following data:

Number of Readings: 327
 Sample Range: 101 to 390 mg/dL

Total within +/- 20 %	
314/327	(96,0%)

Precision

1 lot of Wellion LEONARDO Cholesterol Test Strip5 was used for within-run repeatability study. Venous blood in heparin-tubes was spiked to 3 kinds of concentration. Blood cholesterol readings were recorded for 1 day resulting in 100 data points for each concentration; as shown in the following tables:

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© MED TRUST, Wellion and LEONARDO are registered trademarks.