Disinfectant pad: 70% Alcohol
 Gasa desinfetante: 70% Alcohol
 Gaze desinfetante: 70% Álcool

Vitrex Medical A/S
 Vasekaer 6-8
 DK-2730 Herlev Denmark

Lancet / Lanceta
 Vitrex Medical A/S
 Vasekaer 6-8
 DK-2730 Herlev Denmark

SYMBOLS / SÍMBOLOS / SIMBOLOGIA

IVD In vitro diagnostic device
 Producto sanitario diagnóstico "in vitro"
 Dispositivo médico de diagnóstico In Vitro

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🏭 Legal manufacturer
 Fabricante
 Fabricante

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Registered trademarks.

IRON AND ANAEMIA

Anaemia due to iron depletion is widely held in children and women of all ages but mainly in women who still have their period (at least 20% suffer for iron deficiency). Main signs as paleness, feeling tired, headaches, faster heartbeat, or shortness of breath during exercise appear gradually and could go unnoticed.

Iron deficiency occurs when blood does not contain enough red cells and thus low levels of haemoglobin which is the major protein involved in oxygen transport in whole body. An important component of haemoglobin is iron.

Depletion of iron, which can happen during pregnancy, growth, in case of insufficient iron intake, inadequate absorption or blood loss (period, abnormal bleedings, ulcers, etc.) has tremendous effects on health.

The **IRON-FER TEST** is a rapid immunodiagnostic test for the assessment of ferritin (a protein capable of storing iron in cells) from finger prick whole blood sample. Therefore the Test can be used for the screening of potential iron deficiency.

PRIMA[®]
HOME TEST

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Cut-off: 20 ng/mL

N P

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IRON-FER TEST
TEST DE HIERRO-FER
TESTE DE FERRO-FER

CE
0483

Self-Test for the determination of iron deficiency in blood

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PRECAUTIONS

1) Carefully read the instructions before performing the test. The test is only interpretable if the instructions are carefully respected.

2) Keep out of the reach of children.

3) Do not use after the expiry date printed on the label and on the protective pouch or if the pouch is damaged.

4) Follow strictly the indicated time, blood and diluent quantities.

5) Store between +4 °C and +30 °C. Do not freeze.

6) Do not re-use the cassette and lancet.

7) External use only. DO NOT SWALLOW.

8) This test is exclusively intended to in vitro diagnostic self-testing use.

9) Not recommended for use by individuals taking medicines that thin the blood (anticoagulants) or by the individuals with disorders such as hemophilia.

10) After use, dispose of all components in accordance with local regulations, seek advice from your pharmacist.

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CONTENT

- 1 sealed aluminium pouch containing:

1 plastic pipette

- 2 sterile lancets for blood sampling

- 1 disinfectant pad

1 Ferritin FER Test cassette
 Blue cap
 White cap
 - 1 dropper bottle containing 1 ml of diluent
 Ferritin FER Test

* Not provided necessary material: absorbent cotton.
 * Only open the protective pouch when you are ready to use the test.
 * The desiccant bag should not be used.

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






PROCEDURE

- 1) Wash your hands with soap and rinse with clear warm water.
- 2) Tear the protective pouch (from the notch) and only get out the cassette and the pipette. Dispose of the small desiccant bag.
- 3) Hold the lancet and carefully rotate 360° the protective cap without pulling it. - FIG. A
- 4) Extract and dispose the released cap. - FIG. B
- 5) Press the lancet, on the side from where the cap was extracted, against the fingertip (it is advisable the ring finger side). - FIG. C The tip retracts automatically and safely after use.
- 6) Keeping the hand down massage the end that was pricked to obtain a blood drop. - FIG. D
- 7) Without pressing the pipette bulb, put it in contact with the drop. The blood migrates into the pipette through the capillarity to the line indicated on the pipette. - FIG. E
You may massage again your finger to obtain more blood if the line is not reached. As far as possible, avoid of air bubbles.
- 8) Put the blood collected with the pipette into the sample well of the cassette, by pressing on the pipette bulb. - FIG. F
- 9) Wait for the blood to be totally dispensed in the well. Unscrew the blue cap of the plastic dropper (leave the white cap tightly screwed) and add 5 drops of diluent into the sample well of the cassette. - FIG. G
- 10) Read the result after 10 minutes.

RESULT INTERPRETATION


READ THE RESULT AFTER 10 MINUTES.
DO NOT READ THE RESULT AFTER 15 MINUTES.

The intensity of the line colour does not have any importance for the interpretation of the test result.


NEGATIVE RESULT

Two coloured lines appear in the window under the T (Test) and C (Control) marks. The intensity of the T line may be clearer than the intensity of the C line. This result means that the Ferritin concentration in blood is normal and that there is no potential iron deficiency.


POSITIVE RESULT

Only one coloured line appears under the C (Control) mark. This result means that the ferritin concentration in blood is too low. Reserves are insufficient. You should consult a physician because it may be an iron deficiency.

NON VALID RESULT

No line appears or a coloured line appears under the T (Test) mark without any line under the C (Control) mark. In this case, it is not possible to interpret the test, which must be considered as non valid. It is recommended to repeat the test with a new IRON-FER TEST and a new blood sample.

F.A.Q. - QUESTIONS & ANSWERS

HOW DOES THE IRON-FER TEST WORK? Ferritin is a protein and the primary form of iron stored inside cells. A positive result means that the ferritin concentration in blood is lower than 20 ng/mL and a possible iron deficiency.

WHEN SHOULD THE TEST BE USED? The IRON-FER TEST can be performed in case of symptoms like paleness, feeling tired, headaches, faster heartbeat or shortness of breath during exercise; mainly, if woman, when pregnant or in case of copious periods. The test can be performed anytime of the day but must not be performed in case of disease, acute inflammations or in case of spleen or liver injury. Positive results can be obtained even in case of no iron deficiency situation.

CAN THE RESULT BE INCORRECT? The results are accurate as long as the instructions are carefully respected. Nevertheless, the result can be incorrect if the IRON-FER TEST gets wet before test performing or if the quantity of blood dispensed in the sample well is not sufficient, or if the number of diluent drops are less than 4 or more than 5. The plastic pipette provided in the box allows making sure the collected blood volume is correct.

HOW TO INTERPRET THE TEST IF THE COLOUR AND THE INTENSITY OF THE LINES ARE DIFFERENT? The colour and intensity of the lines have no importance for result interpretation. The lines should only be homogeneous and clearly visible. The test should be considered as negative whatever the colour intensity of the test line is.

WHAT IS THE LINE THAT APPEARS UNDER THE C (CONTROL) MARK FOR? When this line appears, it only means that the test was performing well.

IF I READ THE RESULT AFTER 15 MINUTES, WILL THE RESULT BE RELIABLE? No. The result should be read within 10 minutes after adding the diluent. The result is reliable up to 15 minutes.

WHAT DO I HAVE TO DO IF THE RESULT IS POSITIVE? If the result is positive, it means that the ferritin level is lower than the norm (20 ng/mL) and that you should consult the physician and show him the test. Then, the physician will decide whether additional analysis should be performed.

WHAT DO I HAVE TO DO IF THE RESULT IS NEGATIVE? If the result is negative, it means that the ferritin level is higher than 20 ng/mL and is within the norm. However, if the symptoms persist, it is recommended to consult a physician.

WHAT IS THE ACCURACY OF THE IRON-FER TEST? The IRON-FER TEST is very accurate and has been used for more than 10 years by professionals (hospitals, laboratories, ...) in the field. Evaluation reports show an overall agreement of at least 98% with reference methods.

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